

# PCHi-C Time course Nutlin 3A

---

TUESDAY, 7/20/2021

HCT116 WT cultivated with Nutlin 3A BR2 ([detailed in 210720 Nultin3a timecourse for Hi-C](#)), harvested and fixed for PCHi-C

1. Cells were trypsinized
2. Cells were fixed in suspension with DMEM 10%FBS 2% formaldehyde rotating 10min RT
3. Add ice cold Glycine to 0.125M and rotate 5min RT
4. Incubate 15min on ice with occasional mixing by inversion
5. Centrifuge 800g 10min 4°C
6. Remove supernatant, wash with ice cold 1x PBS
7. Centrifuge 800g 10min 4°C
8. Remove supernatant, flash freeze

TUESDAY, 11/30/2021

Start 1st week Hi-C

Prepare lysis buffer:

- 4690ul ddH2O
- 50ul 1M Tris
- 50ul 20% NP-40/IGEPAL
- 10ul NaCl 5M
- 200ul cOmplete protease inhibitor 25X

1. **Resuspend cells in 1ml lysis buffer**
2. **Incubate 30min on ice** with slight vortex every 5min
3. **Centrifuge them 1000g 10min 4°C**
4. **Remove supernatant**
5. **Resuspend nuclei in 500ul 1.25X NEB2 buffer**
6. **Centrifuge them 1000g 10min 4°C**
7. **Remove supernatant** leaving behind 20ul
8. **Resuspend nuclei in 358ul 1.25X NEB2 Buffer**
9. **Add 11ul SDS 10%** (removes proteins that were not directly crosslinked to the DNA)
10. **Incubate 1h 37°C 950rpm** in a thermomixer
11. **Add 75ul Triton X-100 10%** to quench SDS
12. **Incubate 1h 37°C 950rpm** in a thermomixer
13. **Add 15ul of HindIII** (NEB 100U/ul, R0104T) to digest chromatin
14. **Incubate overnight 37°C 950rpm** in a thermomixer

WEDNESDAY, 12/1/2021

1. **Add an extra 5ul of HindIII** (NEB 100U/ul R0104T) to further digest chromatin
2. **Incubate 1h 37°C 950rpm** in a thermomixer
3. **Prepare the following mastermix (30ul per sample) and add it:**
  - a. 6ul 10x NEB2 Buffer
  - b. 2ul ddH2O
  - c. 1.5ul 10mM dCTP
  - d. 1.5ul 10mM dTTP

- e. 1.5ul 10mM dGTP
- f. 37.5ul 0.4mM biotin-14-dATP (Invitrogen 19524-016)
- g. 10ul Klenow DNA pol I large fragment (NEB 5U/ul M0210L)
4. Mix without creating bubbles and **incubate 75min at 37°C** in a thermomixer. Mix by inversion a couple of times every few minutes.
5. Cool thermomixer to 16°C. **Prepare the following mastermix and add it:**
  - a. 100ul T4 DNA Ligase Buffer 10x (NEB B0202S)
  - b. 5ul BSA 20mg/ml (NEB B9000S)
  - c. 25ul T4 DNA Ligase (Invitrogen 1U/ul 15224-025)
  - d. 346ul ddH2O
6. **Incubate samples 4h 16°C** in a thermomixer. Every now and then mix by inversion.
7. **Incubate samples 30min RT**
8. **Add 10ul RNaseA 10mg/ml and incubate 30min 37°C**
9. **Add 60ul Proteinase K 10mg/ml and incubate 30min 37°C and overnight 65°C 950rpm** in a thermomixer

## THURSDAY, 12/2/2021

1. **Add an extra 30ul of Proteinase K 10mg/ml** to further decrosslink
2. **Incubate 2h 37°C** in a thermomixer
3. **Freeze -20°C** due to time constraints. Phenol-chloroform will be performed in the next days

## THURSDAY, 12/9/2021

1. **Add to the sample 1 volume (545ul) of Phenol:Chloroform**
2. **Vortex** vigorously and **spin 12000g 5min** in a tabletop centrifuge
3. **Transfer the aqueous phase** to a new 2ml tube
4. To precipitate DNA **add the following to the sample** (aprox 545ul):
  - a. 1362.5ul of ice cold EtOH 100%
  - b. 54.5ul sodium acetate 3M pH 5.2
  - c. 2ul Glycoblue
5. **Incubate 1h -80°C**
6. **Spin** sample at **maximum speed 30min 4°C** in a tabletop centrifuge
7. **Discard supernatant** and **add 1.5ml** of freshly made **EtOH 70%**
8. **Spin** sample at **maximum speed 10min RT**
9. **Discard supernatant and air dry pellet**
10. **Resuspend in 83ul TLE**
11. **Quantify** by using Qubit dsDNA HS

Table1



	A	B	C
1	Sample	ng/ul	ng TOTAL
2	HCT116 1h	1370	113710
3	HCT116 4h	1340	111220
4	HCT116 7h	1200	99600
5	HCT116 24h	620	51460

**4h** 0'5ul Dekker // 1ul Roger, 03, Hist

**7h** 0'5ul Dekker // 1ul Roger, 03, Hist

**24h** 1ul Dekker // 2ul Roger, 03, Hist

Tempase Hotstart Pol VWR

95°C 15'

94°C 30" |

60°C 1' | 34 cycles

72°C 1' |

72°C 10'

4°C Hold

[1h] Dekker, Roger, 03, Hist // [4h] Dekker, Roger, 03, Hist // [7h] Dekker, Roger, 03, Hist // [24h] Dekker, Roger, 03, Hist //

[NTC] Dekker, Roger, 03, Hist



Nice. Reamplify 5 cycles 2.5ul PCR product from Dekker

Digest with HindIII, NheI, x2 or nothing

[1h] HindIII, NheI, x2, none // [4h] HindIII, NheI, x2, none // [7h] HindIII, NheI, x2, none // [24h] HindIII, NheI, x2, none



Nice

WEDNESDAY, 12/15/2021

Get 10ug from each sample and top up with ddH2O to 164ul

Table2

	A	B	C
1	Sample	ul library to 10ug	ul ddH2O
2	1h	7.3	156.7
3	4h	7.46	156.54
4	7h	8.33	155.67
5	24h	16.13	147.87

- Add to each 164ul the following mastermix (total volume 200ul)
  - 1ul BSA 20mg/ml
  - 1ul ddH2O
  - 20ul NEB2 10x
  - 4ul dATP
  - 10ul T4 DNA Pol (M0203L)
- Incubate 4h at 20°C no rotation
- Add 4ul EDTA 0.5M and mix
- Top up to 300ul with TLE and perform a Phenol-Chloroform extraction
- Add to the 300ul aqueous phase:
  - 750ul 100% EtOH ice cold
  - 30ul sodium acetate 3M pH 5.2
  - 2ul Glycoblue
- Incubate 1h -80°C
- Centrifuge 4°C 30min max speed
- Discard supernatant
- Resuspend and wash in 1ml EtOH 70%
- Centrifuge RT 10min max speed
- Discard supernatant and resuspend in 130ul

**THURSDAY, 12/16/2021**

1. Transfer to a Covaris 130ul cuvette
2. Sonicate in Covaris M220
  - a. DF 20%
  - b. CPB 200
  - c. PIP 50
  - d. 65sec
3. Transfer to a new 1.5ml low-bind tube
4. Add the following mastermix to each sample:
  - a. 18ul 10x NEB ligation buffer
  - b. 18ul dNTPs 2.5mM each
  - c. 6.5ul T4 DNA Pol (NEB M0203)
  - d. 6.5ul T4 PNK (NEB M0201)
  - e. 1.3ul Klenow (NEB M0210)
5. Incubate 20°C 30min
6. Add 130ul SPRI beads (1vol) and perform standard DNA cleanup
7. Resuspend in 35.7ul TLE

**FRIDAY, 12/17/2021**